

Dry Roll-to-Roll Forming Laboratory Machine



Product Overview

- The equipment is used to realize the film formation and thinning of dry electrode powder and particles

Product Features

- The equipment is sophisticated, and the equipment placement is not limited by the venue
- The operation is convenient and meets the operational needs of the glove box
- The double rollers are driven independently to realize the flexible operation of the equipment

Product Parameters

Mechanical speed	MAX:10m/min
Operation speed	MAX:5m/mim
Roller width	200mm
Effective roller width	120mm
Operating temperature	≤200°C
Rolling pressure	1T on one side, 2T in total,display accuracy 0.1%
Calendering thickness accuracy	±10μm